

张家港康百隆机械有限公司
Zhangjiagang Combilong Machinery Co.,Ltd

Proforma Invoice

InvoiceNO.: 2020072709

Date: 2020/07/27

Supplier: Zhangjiagang Combilong Machinery Co.,Ltd

Contact Person: Melissa Miao

E-mail: melissa@combimedical.com

Buyer Company:

Contact Person:




Phone:

FROM: CHINA







TO: Lima

No.	Name	Photo	Qty	Unit Price (USD)	Total Price (USD)	Remark
1	Isolator Stretcher		6 sets	2,850	17,100	
<i>Optional</i>	<i>Hepa Filter (Optional)</i>		30 set (4 pcs)	80	2,400	<i>disposable</i>
2	Feight				1,807	
3	Tax				3,510	
Total Amount	US DOLLARS TWENTY-FOUR THOUSAND,EIGHT HUNDRED AND SEVENTEEN ONLY				24,817	
Remarks	1. This is CIF Shanghai price. 2. Payment term: 100% deposit. 3. Date of shipment: Arrange shipment within 3 working days once receiving payment of deposit. 4. Packing: standard carton packing.					
BANK ACCOUNT	Bank Details: Bank of China, Zhangjiagang sub-branch Swift No.: BKCHCNBJ95L Account No.: 468960006190 Beneficiary: Zhangjiagang Combilong Machinery Co., Ltd. Bank Address: Zhangjiagang City, Jiangsu Province, China					


张家港康百隆机械有限公司
 Zhangjiagang Combilong Machinery Co.,Ltd

<u>Product description</u>	Isolator Stretcher for Infected Patient
<u>Product photo</u>	 <p style="text-align: center;">Oxygen supply connector</p>
<u>Operating Principle</u>	<p>With negative pressure system, inside the chamber it will remain a stable negative-pressure environment and prevent the polluted air inside from leaking out so that we can keep the medical workers safe.</p> <p>In addition, a suction blower will operate to pump the inside air circulatingly after being cleaned by the hepa filters at the outlet to keep the inside air fresh.</p>
<u>Overall Structure</u>	
<u>Photos</u>	

张家港康百隆机械有限公司
Zhangjiagang Combilong Machinery Co.,Ltd

		
		
		
<u>Specs</u> Dimintions:		
Height	430mm	
Width	570mm	
Length	1950mm	
Negative Pressure of Blower	> 20Pa	
Air Flow	>100L/min	

张家港康百隆机械有限公司
Zhangjiagang Combilong Machinery Co.,Ltd

Volts	Battery: 36V.4Ah	
Filtration efficiency	0.3micra > 99.9995 %	
Consum	30 Va	
Battery autonomy	> 5 hours	
Weight	17kg	
Main accessory	<div> <div>Filter</div> <div>4pc</div> </div> <div> <div>Glove</div> <div>1 pair</div> </div> <div> <div>Airway tube</div> <div>1 set</div> </div> <div> <div>Battery and recharger</div> <div>1 set</div> </div> <div> <div>Motor case and work fan</div> <div>1set</div> </div>	
Note	The filter is disposable, 4 pcs/set, they should be changed after one patient used. And it shall be used within 12 hours after unpacking.	